

Anti-Non-Phospho (Active) B-Catenin (S33/S37/T41) Monoclonal Antibody For Wb, Ihc-P, Elisa - P35222

Item Number: CM0006557



Introduction

KO validated rabbit monoclonal antibody specific for non-phosphorylated (active) β -catenin at S33, S37, T41. Suitable for WB, IHC-P, ELISA; detects human, mouse, rat β -catenin (85 kDa). Ideal for Wnt signaling, EMT, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Non-phospho (Active) β -Catenin -S33/S37/T41 Rabbit mAb
Remarks/Alias	EV7; CTNNB; MRD19; NEDSDV; armadillo; Non-phospho (Active) β -Catenin -S33/S37/T41
Species	Human
GeneID (Human)	1499
GeneID	1499
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Non-phospho (Active) β -Catenin -S33/S37/T41 (NP_001895.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag