


Anti-Nono/P54Nrb Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15233

Item Number: CM0027888



Introduction

Rabbit polyclonal antibody targeting human NONO/p54nrb (Q15233), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for studying RNA splicing, DNA repair, paraspeckle biology, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	NONO/p54nrb Rabbit pAb
Remarks/Alias	P54; NMT55; NRB54; MRXS34; P54NRB; PPP1R114; NONO/p54nrb
Species	Human
Gene ID (Human)	4841
Gene ID	4841
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NONO/p54nrb (NP_001138881.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15233
Protein Weight	54kDa
Shipping	Ice bag