


Anti-Nono/P54Nrb Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15233

Item Number: CM0008374



Introduction

NONO/p54nrb rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat NONO/p54nrb (Q15233). Ideal for splicing, paraspeckle, and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	NONO/p54nrb Rabbit mAb
Remarks/Alias	P54; NMT55; NRB54; MRXS34; P54NRB; PPP1R114; NONO/p54nrb
Species	Human
GeneID (Human)	4841
GeneID	4841
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15233
Protein Weight	54 kDa
Shipping	Ice bag