


Anti-Nop10 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Npe3

Item Number: CM0010132



Introduction

Rabbit monoclonal anti-NOP10 antibody validated for WB, IHC-P, IP, and ELISA. Cross-reacts with human, mouse, and rat NOP10. Ideal for studying ribosome biogenesis, telomere maintenance, and pseudouridylation.

[Learn More](#)

Parameter	Value
Product Name	NOP10 Rabbit mAb
Remarks/Alias	DKCB1; NOLA3; NOP10P; NOP10
Species	Human
Gene ID (Human)	55505
Gene ID	55505
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-64 of human NOP10 (Q9NPE3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPE3
Protein Weight	8kDa
Shipping	Ice bag