


Anti-Nop14 Polyclonal Antibody For Wb, Elisa - P78316

Item Number: CM0017168



Introduction

NOP14 rabbit polyclonal antibody validated for WB and ELISA. Targets human nucleolar protein 14, crucial for pre-18S rRNA processing and 40S pre-ribosomal export. Ideal for ribosome biogenesis research.

[Learn More](#)

Parameter	Value
Product Name	NOP14 Rabbit pAb
Remarks/Alias	UTP2; NOL14; C4orf9; RES425; RES4-25; NOP14
Species	Human
GeneID (Human)	8602
GeneID	8602
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78316
Protein Weight	98kDa
Shipping	Ice bag