


Anti-Nop16 Polyclonal Antibody For Wb, Elisa - Q9Y3C1

Item Number: CM0014961



Introduction

NOP16 Rabbit Polyclonal Antibody, validated for WB and ELISA. Recognizes Human, Mouse, Rat NOP16. Recombinant immunogen. Ideal for nucleolar protein research. 21 kDa.

[Learn More](#)

Parameter	Value
Product Name	NOP16 Rabbit pAb
Remarks/Alias	HSPC111; HSPC185; NOP16
Species	Human
GeneID (Human)	51491
GeneID	51491
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human NOP16 (NP_057475.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3C1
Protein Weight	21kDa
Shipping	Ice bag