


Anti-Nop2 Polyclonal Antibody For Wb, Elisa - P46087

Item Number: CM0023336



Introduction

NOP2 Rabbit pAb, a polyclonal antibody targeting human NOP2, validated for WB and ELISA with human/mouse cross-reactivity. Ideal for rRNA processing and ribosome biogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	NOP2 Rabbit pAb
Remarks/Alias	NOL1; p120; NSUN1; NOP120; NOP2
Species	Human
GeneID (Human)	4839
GeneID	4839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 355-590 of human NOP2 (NP_001245237.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46087
Protein Weight	89kDa
Shipping	Ice bag