


Anti-Nop56 Polyclonal Antibody For Wb, Elisa - 000567

Item Number: CM0021647



Introduction

Rabbit polyclonal antibody against human NOP56 (NOL5A, SCA36). Validated for Western blot (WB) and ELISA. Ideal for studies on nucleolar protein function and ribosome biogenesis.

[Learn More](#)

Parameter	Description
Product Name	NOP56 Rabbit pAb
Remarks/Alias	NOL5A; SCA36; NOP56
Species	Human
GeneID (Human)	10528
GeneID	10528
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	000567
Protein Weight	66kDa
Shipping	Ice bag