

Anti-Nop58 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y2X3

Item Number: CM0017155



Introduction

NOP58 Rabbit pAb is a polyclonal antibody raised against a recombinant fusion protein, validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for ribosomal biogenesis and nucleolar research.

[Learn More](#)

Parameter	Value
Product Name	NOP58 Rabbit pAb
Remarks/Alias	NOP5, HSPC120, NOP5/NOP58, NOP58
Species	Human
GeneID	51602
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-529 of human NOP58 (NP_057018.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2X3
Protein Weight	60kDa
Shipping	Ice bag