

Anti-Nos1 Polyclonal Antibody For Wb & Elisa - P29475

Item Number: CM0014879



Introduction

High-quality rabbit polyclonal antibody against human NOS1 (neuronal nitric oxide synthase) for Western blot and ELISA applications. Cross-reacts with rat. Produces nitric oxide, a key neurotransmitter and signaling molecule. Ideal for neuroscience and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	NOS1 Rabbit pAb
Remarks/Alias	NOS; bNOS; nNOS; IHP51; N-NOS; NC-NOS; NOS1
Species	Human
Gene ID (Human)	4842
Gene ID	4842
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P29475
Protein Weight	161kDa
Shipping	Ice bag