

Anti-Nos1 Polyclonal Antibody For Wb, If/Icc, Elisa - P29475

Item Number: CM0016401



Introduction

Rabbit polyclonal antibody targeting human NOS1 (neuronal nitric oxide synthase), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat. For research on nitric oxide signaling, synaptic transmission, and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	NOS1 Rabbit pAb
Remarks/Alias	NOS; bNOS; nNOS; IHPS1; N-NOS; NC-NOS; NOS1
Species	Human
GeneID (Human)	4842
GeneID	4842
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human NOS1 (NP_000611.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P29475
Protein Weight	161kDa
Shipping	Ice bag