


Anti-Nosip Polyclonal Antibody For Wb, Elisa - Q9Y314

Item Number: CM0017167



Introduction

High-quality rabbit polyclonal antibody against human NOSIP (eNOS-interacting protein), validated for WB and ELISA. Ideal for studies on nitric oxide regulation and developmental biology. Cross-reacts with human. Protein weight: 33 kDa.

[Learn More](#)

Parameter	Value
Product Name	NOSIP Rabbit pAb
Remarks/Alias	CGI-25; NOSIP
Species	Human
GeneID (Human)	51070
GeneID	51070
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human NOSIP (NP_057037.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y314
Protein Weight	33kDa
Shipping	Ice bag