


Anti-Nostrin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ivi9

Item Number: CM0023394



Introduction

CamelBio offers high-quality NOSTRIN Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Specific to human NOSTRIN with cross-reactivity to mouse and rat. Ideal for eNOS trafficking and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	NOSTRIN Rabbit pAb
Remarks/Alias	DaIP2; NOSTRIN
Species	Human
GeneID (human)	115677
GeneID	115677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-430 of human NOSTRIN (NP_001034813.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8Ivi9
Protein Weight	58kDa
Shipping	Ice bag