

Anti-Notch1 Monoclonal Antibody For Wb, Ip, Elisa - P46531

Item Number: CM0003362



Introduction

Rabbit monoclonal antibody against human Notch1 (P46531) validated for WB, IP, ELISA. Cross-reacts with mouse. Ideal for studying Notch signaling, cell differentiation, and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	Notch1 Rabbit mAb
Remarks/Alias	hN1; AOS5; TAN1; AOVD1; Notch1
Species	Human
GeneID (Human)	4851
GeneID	4851
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 2456-2555 of human Notch1 (P46531).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46531
Protein Weight	273kDa
Shipping	Ice bag