

## Anti-Notch1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P46531

Item Number: CM0029394



### Introduction

Rabbit polyclonal antibody against human Notch1 (neurogenic locus notch homolog protein 1), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying Notch signaling and cell-fate determination.

[Learn More](#)

| Parameter        | Value                          |
|------------------|--------------------------------|
| Product Name     | Notch1 Rabbit pAb              |
| Remarks/Alias    | hN1; AOS5; TAN1; AOVD1; Notch1 |
| Species          | Human                          |
| GeneID (Human)   | 4851                           |
| GeneID           | 4851                           |
| Immunogen        | Recombinant protein            |
| Source           | Rabbit                         |
| Category         | Polyclonal Antibodies          |
| Application      | WB, IHC-P, IF/ICC, ELISA       |
| Cross Reactivity | Human, Mouse, Rat              |
| SWISS            | P46531                         |
| Protein Weight   | 273kDa                         |
| Shipping         | Ice bag                        |