

Anti-Notch1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P46531

Item Number: CM0019838



Introduction

Rabbit polyclonal antibody against Notch1 (aa 2428-2555) validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat Notch1. Suitable for studying Notch signaling in development and disease.

[Learn More](#)

Parameter	Value
Product Name	Notch1 Rabbit pAb
Remarks/Alias	hN1; AOS5; TAN1; AOVD1; Notch1
Species	Human
GeneID (Human)	4851
GeneID	4851
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2428-2555 of human Notch1 (NP_060087.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46531
Protein Weight	273kDa
Shipping	Ice bag