

Anti-Notch2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q04721

Item Number: CM0014118



Introduction

Rabbit polyclonal antibody against human NOTCH2, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for studying Notch signaling, cell fate determination, and osteoclast differentiation.

[Learn More](#)

Parameter	Value
Product Name	NOTCH2 Rabbit pAb
Remarks/Alias	hN2; AGS2; HJCYS; NOTCH2
Species	Human
GeneID (Human)	4853
GeneID	4853
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q04721
Protein Weight	265kDa
Shipping	Ice bag