

Anti-Notch2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q04721

Item Number: CM0023175



Introduction

Anti-NOTCH2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Target: ~265kDa NOTCH2 receptor (Q04721), involved in cell-fate determination, bone remodeling, and cancer.

[Learn More](#)

Parameter	Specification
Product Name	NOTCH2 Rabbit pAb
Remarks/Alias	NOTCH2; AGS2; HJCYS; hN2; notch 2
Species	Human
Gene ID	4853
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1697-1819 of human NOTCH2 (NP_077719.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04721
Protein Weight	265kDa
Shipping	Ice bag