

Anti-Notch3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Um47

Item Number: CM0017423



Introduction

NOTCH3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying Notch signaling, CADASIL, and cell-fate determination.

[Learn More](#)

Parameter	Value
Product Name	NOTCH3 Rabbit pAb
Remarks/Alias	IMF2, LMNS, CASIL, CADASIL, CADASIL1, NOTCH3
Species	Human
Gene ID (Human)	4854
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 751-850 of human NOTCH3 (NP_000426.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UM47
Protein Weight	244kDa
Shipping	Ice bag