

Anti-Notch3 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Um47

Item Number: CM0016182



Introduction

Polyclonal rabbit antibody targeting human NOTCH3 (Q9UM47). Suitable for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for Notch signaling, developmental biology, and CADASIL research.

[Learn More](#)

Parameter	Value
Product Name	NOTCH3 Rabbit pAb
Remarks/Alias	IMF2; LMNS; CASIL; CADASIL; CADASIL1; NOTCH3
Species	Human
Gene ID (Human)	4854
Gene ID	4854
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2200-2321 of human NOTCH3 (NP_000426.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UM47
Protein Weight	244kDa
Shipping	Ice bag