


Anti-Notch4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99466

Item Number: CM0025454



Introduction

Rabbit polyclonal antibody against human NOTCH4 (INT3), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for Notch signaling and developmental biology research. SWISS: Q99466.

[Learn More](#)

Parameter	Value
Product Name	NOTCH4 Rabbit pAb
Remarks/Alias	INT3; NOTCH4
Species	Human
Gene ID (Human)	4855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1824-2003 of human NOTCH4 (NP_004548.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99466
Protein Weight	210kDa
Shipping	Ice bag