

Anti-Nov (Ccn3/Igfbp9) Polyclonal Antibody For Ihc-P And Elisa - P48745

Item Number: CM0027791

Introduction

Rabbit polyclonal antibody targeting human NOV (CCN3/IGFBP9) for IHC-P and ELISA applications. Recognizes a 39 kDa protein involved in cell proliferation, adhesion, angiogenesis, and bone regeneration.

[Learn More](#)

Parameter	Value
Product Name	NOV Rabbit pAb
Remarks/Alias	NOV; NOVh; IBP-9; IGFBP9; IGFBP-9
Species	Human
Gene ID (Human)	4856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-357 of human NOV (NP_002505.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P48745
Protein Weight	39kDa
Shipping	Ice bag