

Anti-Nox1 Polyclonal Antibody For Wb, Elisa - Q9Y5S8

Item Number: CM0030353

Introduction

Rabbit polyclonal antibody targeting human NADPH oxidase 1 (NOX1). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Useful for superoxide signaling, cancer biology, and colon tissue studies.

[Learn More](#)

Parameter	Value
Product Name	NOX1 Rabbit pAb
Remarks/Alias	MOX1; NOH1; NOH-1; GP91-2; NOH-1L; NOX1
Species	Human
Gene ID (Human)	27035
Gene ID	27035
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 241-340 of human NOX1 (NP_008983.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5S8
Protein Weight	65kDa
Shipping	Ice bag