

Anti-Nox1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5S8

Item Number: CM0015750

Introduction

Rabbit polyclonal antibody against human NOX1 (NADPH oxidase 1), validated for WB, IHC-P, ELISA. Cross-reacts with human and mouse. Suitable for studying superoxide generation and colon carcinoma.

[Learn More](#)

Parameter	Value
Product Name	NOX1 Rabbit pAb
Remarks/Alias	MOX1; NOH1; NOH-1; GP91-2; NOH-1L; NOX1
Species	Human
Gene ID (Human)	27035
Gene ID	27035
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 418-564 of human NOX1 (NP_008983.2?).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5S8
Protein Weight	65kDa
Shipping	Ice bag