

Anti-Nox1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y5S8

Item Number: CM0010204

Introduction

Anti-NOX1 Rabbit Monoclonal Antibody (validated for WB and ELISA). Detects human NOX1, a 65 kDa NADPH oxidase. Suitable for colon, uterus, prostate research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] NOX1 Rabbit mAb
Remarks/Alias	MOX1; NOH1; NOH-1; GP91-2; NOH-1L
Species	Human
Gene ID (Human)	27035
Immunogen	Synthetic peptide corresponding to amino acids 241-340 of human NOX1 (NP_008983.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y5S8
Protein Weight	65kDa
Shipping	Ice bag