

Anti-Nox1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5S8

Item Number: CM0014321

Introduction

Anti-NOX1 rabbit polyclonal antibody designed for WB, IF/ICC, and ELISA applications. Recognizes human NADPH oxidase 1 (NOX1), a membrane-bound superoxide generator linked to colon and prostate cancer. Ideal for diagnostic research and assay development.

[Learn More](#)

Parameter	Value
Product Name	NOX1 Rabbit pAb
Remarks/Alias	MOX1; NOH1; NOH-1; GP91-2; NOH-1L; NOX1
Species	Human
Gene ID (Human)	27035
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NOX1 (NP_008983.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y5S8
Protein Weight	65 kDa
Shipping	Ice bag