

Anti-Nox2/Gp91Phox Polyclonal Antibody For Wb, Elisa - P04839

Item Number: CM0016787



Introduction

Rabbit polyclonal antibody against human NOX2/gp91phox (CYBB). Validated for WB and ELISA; cross-reacts with human, mouse, rat. Suitable for phagocyte NADPH oxidase and superoxide studies, including CGD research.

[Learn More](#)

Parameter	Value
Product Name	NOX2/gp91phox Rabbit pAb
Remarks/Alias	CGD; CGDX; NOX2; IMD34; AMCBX2; GP91-1; GP91PHOX; p91-PHOX; GP91-PHOX; NOX2/gp91phox
Species	Human
GeneID (human)	1536
GeneID	1536
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 101-200 of human NOX2/gp91phox (NP_000388.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04839
Protein Weight	65kDa
Shipping	Ice bag