

Anti-Nox2/Gp91Phox Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04839

Item Number: CM0022572



Introduction

Rabbit polyclonal antibody against human NOX2/gp91phox, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Research use for studying NADPH oxidase complex and superoxide generation.

[Learn More](#)

Parameter	Value
Product Name	NOX2/gp91phox Rabbit pAb
Remarks/Alias	CGD; CGDX; NOX2; IMD34; AMCBX2; GP91-1; GP91PHOX; p91-PHOX; GP91-PHOX; NOX2/gp91phox
Species	Human
GeneID (Human)	1536
GeneID	1536
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 385-570 of human NOX2/gp91phox (NP_000388.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04839
Protein Weight	65kDa
Shipping	Ice bag