

Anti-Nox2/Gp91Phox Rabbit Monoclonal Antibody For Wb And Elisa - P04839

Item Number: CM0002209



Introduction

High-quality anti-NOX2/gp91phox rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat NOX2. Ideal for studying NADPH oxidase function in phagocytes.

[Learn More](#)

Parameter	Value
Product Name	NOX2/gp91phox Rabbit mAb
Remarks/Alias	CGD; CGDX; NOX2; IMD34; AMCBX2; GP91-1; GP91PHOX; p91-PHOX; GP91-PHOX; NOX2/gp91phox
Species	Human
GeneID (Human)	1536
GeneID	1536
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04839
Protein Weight	65kDa
Shipping	Ice bag