

Anti-Noxa2/P67Phox Polyclonal Antibody For Wb, If-P, Elisa - P19878

Item Number: CM0014730



Introduction

Rabbit polyclonal antibody targeting human NOXA2/p67phox (NCF2) for WB, IF-P, ELISA. Cross-reacts with mouse. Recombinant immunogen; suitable for NADPH oxidase and superoxide production research.

[Learn More](#)

Parameter	Value
Product Name	NOXA2/p67phox Rabbit pAb
Remarks/Alias	NCF-2; NOXA2; P67PHOX; P67-PHOX; NOXA2/p67phox
Species	Human
GeneID (Human)	4688
GeneID	4688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 227-526 of human NOXA2/p67phox (NP_000424.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19878
Protein Weight	60kDa
Shipping	Ice bag