

Anti-Noxa2/P67Phox Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P19878

Item Number: CM0012020



Introduction

High-quality rabbit monoclonal antibody targeting human NOXA2/p67phox (NCF2), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for phagocyte oxidase research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	NOXA2/p67phox Rabbit mAb
Remarks/Alias	NCF-2, NOXA2, P67PHOX, P67-PHOX, NOXA2/p67phox
Species	Human
GeneID (Human)	4688
GeneID	4688
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human NOXA2/p67phox (P19878).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P19878
Protein Weight	60kDa
Shipping	Ice bag