


Anti-Noxo1 Polyclonal Antibody For Wb, Elisa - Q8Nfa2

Item Number: CM0013543



Introduction

Rabbit polyclonal antibody against human NOXO1 (NADPH oxidase organizer 1). Validated for WB and ELISA with cross-reactivity to mouse. Ideal for NOX complex and superoxide generation studies.

[Learn More](#)

Parameter	Value
Product Name	NOXO1 Rabbit pAb
Remarks/Alias	SNX28; P41NOX; P41NOXA; P41NOXB; P41NOXC; SH3PXD5; NOXO1
Species	Human
GeneID	124056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 121-370 of human NOXO1 (NP_653204.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8NFA2
Protein Weight	41kDa
Shipping	Ice bag