

Anti-Npas2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99743

Item Number: CM0019394

Introduction

NPAS2 Rabbit pAb (Q99743) validated for WB, IF/ICC, and ELISA. Detects human NPAS2 with cross-reactivity to mouse and rat. Ideal for circadian clock and transcriptional regulation research.

[Learn More](#)

Parameter	Specification
Product Name	NPAS2 Rabbit pAb
Remarks/Alias	MOP4; PASD4; bHLHe9; NPAS2
Species	Human
GeneID (Human)	4862
GeneID	4862
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 45-240 of human NPAS2 (NP_002509.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99743
Protein Weight	92kDa
Shipping	Ice bag