

Anti-Npas3 Polyclonal Antibody For Wb And Elisa - Q8lxf0

Item Number: CM0019573

Introduction

Rabbit polyclonal antibody targeting human NPAS3, validated for WB and ELISA applications. Detects human and rat NPAS3 (~101 kDa). Ideal for neurogenesis and schizophrenia research.

[Learn More](#)

Parameter	Value
Product Name	NPAS3 Rabbit pAb
Remarks/Alias	MOP6; PASD6; bHLHe12; NPAS3
Species	Human
GeneID (human)	64067
GeneID	64067
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human NPAS3 (NP_001158221.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IXF0
Protein Weight	101kDa
Shipping	Ice bag