

Anti-Npb Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Ng41

Item Number: CM0020207

Introduction

Rabbit polyclonal antibody against human Neuropeptide B, validated for WB, IF/ICC, ELISA. Reacts with human and mouse. Suitable for neuroscience and cancer research applications, including pain and feeding studies.

[Learn More](#)

Parameter	Value
Product Name	NPB Rabbit pAb
Remarks/Alias	L7; PPL7; PPNPB; NPB
Species	Human
GeneID (Human)	256933
GeneID	256933
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-125 of human NPB (NP_683694.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NG41
Protein Weight	13kDa
Shipping	Ice bag