

Anti-Npc2 (He1, Eddm1) Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P61916

Item Number: CM0028328



Introduction

Rabbit polyclonal antibody targeting human NPC2 (HE1, EDDM1), validated for WB, IF-P, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for studying cholesterol transport and Niemann-Pick disease type C2.

[Learn More](#)

Parameter	Value
Product Name	NPC2 Rabbit pAb
Remarks/Alias	HE1; EDDM1; NPC2
Species	Human
GeneID (Human)	10577
GeneID	10577
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-151 of human NPC2 (NP_006423.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61916
Protein Weight	17kDa
Shipping	Ice bag