

Anti-Npc2 Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - P61916

Item Number: CM0007627



Introduction

Recombinant rabbit monoclonal antibody targeting human NPC2, validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for cholesterol transport and Niemann-Pick disease type C studies.

[Learn More](#)

Parameter	Value
Product Name	Niemann-Pick type C2 NPC2 Rabbit mAb
Remarks/Alias	HE1; EDDM1; Niemann-Pick type C2 NPC2
Species	Human
Gene ID (Human)	10577
Gene ID	10577
Immunogen	Recombinant protein of human Niemann-Pick type C2 NPC2
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61916
Protein Weight	17kDa
Shipping	Ice bag