


Anti-Npff Polyclonal Antibody For Wb, Elisa - O15130

Item Number: CM0025786



Introduction

High-quality rabbit polyclonal antibody against human NPFF, validated for Western blot and ELISA. Cross-reacts with rat, suitable for neuropeptide function and pain modulation studies.

[Learn More](#)

Parameter	Value
Product Name	NPFF Rabbit pAb
Remarks/Alias	FMRFAL; NPFF
Species	Human
GeneID (human)	8620
GeneID	8620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-113 of human NPFF (NP_003708.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O15130
Protein Weight	12kDa
Shipping	Ice bag