


Anti-Nphp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15259

Item Number: CM0029171



Introduction

High-quality rabbit polyclonal anti-NPHP1 antibody validated for WB, IHC-P, IF/ICC, and ELISA; detects human, mouse, and rat nephrocystin-1; ideal for ciliopathy and kidney disease research.

[Learn More](#)

Parameter	Value
Product Name	NPHP1 Rabbit pAb
Remarks/Alias	NPH1; JBTS4; SLSN1; NPHP1
Species	Human
Gene ID (Human)	4867
Gene ID	4867
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15259
Protein Weight	83kDa
Shipping	Ice bag