

Anti-Nphp4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75161

Item Number: CM0030783

Introduction

High-quality rabbit polyclonal antibody targeting human NPHP4 (nephrocystin-4) for Western blot, immunofluorescence, and ELISA. Reacts with human, mouse, and rat samples, supporting ciliopathy and renal disease research.

[Learn More](#)

Parameter	Value
Product Name	NPHP4 Rabbit pAb
Remarks/Alias	POC10; SLSN4; NPHP4
Species	Human
GeneID (human)	261734
GeneID	261734
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1247-1426 of human NPHP4 (NP_055917.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75161
Protein Weight	158kDa
Shipping	Ice bag