

Anti-Nphs2 Polyclonal Antibody For Wb, Elisa - Q9Np85

Item Number: CM0020171

Introduction

High-quality rabbit polyclonal antibody against human NPHS2 (podocin). Validated for WB and ELISA, cross-reacts with human, mouse, rat. Essential tool for glomerular permeability and nephrotic syndrome research.

[Learn More](#)

Parameter	Value
Product Name	NPHS2 Rabbit pAb
Remarks/Alias	PDCN; SRN1; NPHS2
Species	Human
GeneID (Human)	7827
GeneID	7827
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-388 of human NPHS2 (NP_055440.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NP85
Protein Weight	42kDa
Shipping	Ice bag