

Anti-Nphs2 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Np85

Item Number: CM0006699

Introduction

Recombinant rabbit monoclonal antibody targeting human podocin (NPHS2) for Western blot (WB) and ELISA. Expressed in kidney podocytes, podocin regulates glomerular filtration. Ideal for nephrotic syndrome research and diagnostics.

[Learn More](#)

Parameter	Value
Product Name	NPHS2 Rabbit mAb
Remarks/Alias	PDCN; SRN1; NPHS2
Species	Human
GeneID (Human)	7827
GeneID	7827
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 284-383 of human NPHS2 (Q9NP85).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NP85
Protein Weight	42kDa
Shipping	Ice bag