


Anti-Npl Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bxd5

Item Number: CM0025706



Introduction

Anti-NPL rabbit polyclonal antibody for WB and ELISA applications. Detects human N-acetylneuraminase lyase, with cross-reactivity to mouse and rat. Recombinant protein immunogen. Used in sialic acid metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	NPL Rabbit pAb
Remarks/Alias	NAL; C112; NPL1; C1orf13; NPL
Species	Human
Gene ID (Human)	80896
Gene ID	80896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human NPL (NP_110396.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9BXD5
Protein Weight	35kDa
Shipping	Ice bag