


Anti-Nploc4 Rabbit Polyclonal Antibody For Wb - Q8Tat6

Item Number: CM0028003



Introduction

NPLOC4 Rabbit Polyclonal Antibody (WB, IF/ICC, ELISA) targets human NPLOC4, a key component of the ubiquitin-dependent protein degradation pathway. Validated in Human, Mouse, and Rat samples.

[Learn More](#)

Parameter	Details
Product Name	NPLOC4 Rabbit pAb
Remarks/Alias	NPL4; NPLOC4
Species	Human
GeneID (human)	55666
GeneID	55666
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TAT6
Protein Weight	68kDa
Shipping	Ice bag