

Anti-Npm3 Polyclonal Antibody For Wb, Elisa - 075607

Item Number: CM0020832

Introduction

Rabbit polyclonal antibody targeting human NPM3 (Nucleoplasmin-3), validated for Western blot and ELISA applications. NPM3 plays a role in chromatin remodeling and ribosome biogenesis, making this antibody useful for nuclear protein studies.

[Learn More](#)

Parameter	Value
Product Name	NPM3 Rabbit pAb
Remarks/Alias	PORMIN; TMEM123; NPM3
Species	Human
Gene ID (Human)	10360
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-140 of human NPM3 (NP_008924.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75607
Protein Weight	19kDa
Shipping	Ice bag