


Anti-Npr3 (Anp-C) Rabbit Polyclonal Antibody For Wb, Elisa - P17342

Item Number: CM0025605



Introduction

Rabbit polyclonal antibody against human NPR3 (ANP-C clearance receptor). Validated for WB and ELISA, reacts with human, mouse, rat. Ideal for natriuretic peptide signaling research.

[Learn More](#)

Parameter	Value
Product Name	NPR3 Rabbit pAb
Remarks/Alias	NPRC; ANP-C; ANPRC; BOMOS; NPR-C; ANPR-C; GUCY2B; C5orf23; NPR3
Species	Human
GeneID (Human)	4883
GeneID	4883
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17342
Protein Weight	60kDa
Shipping	Ice bag