


Anti-Nprl2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wtw4

Item Number: CM0031045



Introduction

High-quality rabbit polyclonal antibody targeting human NPRL2 for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Essential for GATOR1 complex and mTORC1 pathway research. Predicted molecular weight 44 kDa.

[Learn More](#)

Parameter	Value
Product Name	NPRL2 Rabbit pAb
Remarks/Alias	NPR2; NPR2L; TUSC4; FFEVF2; NPRL2
Species	Human
GeneID (human)	10641
GeneID	10641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human NPRL2 (NP_006536.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8Wtw4
Protein Weight	44kDa
Shipping	Ice bag