

## Anti-Npsr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q6W5P4

Item Number: CM0019229

### Introduction

High-quality rabbit polyclonal antibody targeting human NPSR1 (neuropeptide S receptor). Validated for WB and ELISA, specific for human, and raised against a C-terminal peptide. Supports GPCR signaling and asthma research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | NPSR1 Rabbit pAb                                                                                                          |
| Remarks/Alias    | GPRA; NPSR; VRR1; ASRT2; PGR14; GPR154; NPSR1                                                                             |
| Species          | Human                                                                                                                     |
| GeneID (Human)   | 387129                                                                                                                    |
| GeneID           | 387129                                                                                                                    |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human NPSR1 (NP_001287863.1). |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, ELISA                                                                                                                 |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | Q6W5P4                                                                                                                    |
| Protein Weight   | 43kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |