

Anti-Nptx2 Rabbit Polyclonal Antibody For Wb And Elisa - P47972

Item Number: CM0014069

Introduction

Rabbit polyclonal antibody targeting human neuronal pentraxin-2 (NPTX2), suitable for Western blot and ELISA. Supports research into synaptic plasticity and neurobiology.

[Learn More](#)

Parameter	Value
Product Name	NPTX2 Rabbit pAb
Remarks/Alias	NP2, NARP, NP-II, NPTX2
Species	Human
Gene ID (Human)	4885
Gene ID	4885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-240 of human NPTX2 (NP_002514.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P47972
Protein Weight	47kDa
Shipping	Ice bag