

Anti-Npy Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01303

Item Number: CM0026511



Introduction

Rabbit polyclonal antibody against human NPY (Pro-neuropeptide Y), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for neuropeptide Y signaling and feeding behavior research.

[Learn More](#)

Parameter	Value
Product Name	NPY Rabbit pAb
Remarks/Alias	PYY4; NPY
Species	Human
GeneID (Human)	4852
GeneID	4852
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-97 of NPY (NP_000896.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01303
Protein Weight	11kDa
Shipping	Ice bag